

The Challenge of Detecting Prodromal Alzheimer's Disease for Clinical Trials: Enriching Recruitment Through Identification of Episodic Memory Deficits at Pre-Screening

Rosemary A. Abbott¹, Francesca Cormack¹, Kenton Zavitz¹, Babette L.R. Reijs², Inez H.G.B. Ramakers², Charlotte E. Teunissen³, Marleen Koel-Simmelink³, Magda Tsolaki⁴, Lars-Olof Wahlund⁵, Gunhild Waldemar⁶, Lucrezia Hausner⁷, Rik Vandenberghe⁸, Peter Johannsen⁶, Hugo Vanderstichele⁹, Frans Verhey², Andrew Blackwell¹, Jenny H. Barnett¹ and Pieter Jelle Visser^{2,3}

¹Cambridge Cognition, Cambridge, UK, ²Maastricht University NL, ³VU University Medical Centre, NL, ⁴Aristotle University of Thessaloniki, Greece ⁵Karolinska Institute, Sweden, ⁶Copenhagen University Hospital, Denmark, ⁷Heidelberg University, Germany, ⁸University Hospital, Leuven, Belgium, ⁹ADx NeuroSciences, Belgium.

Background

- Alzheimer's disease (AD) trials require participants with a specific biomarker and cognitive profile, targeting particular disease severity (e.g. prodromal vs MCI vs AD).
- Finding participants for clinical trials in the 'pre-clinical' or 'prodromal' stage of Alzheimer's disease prior to the onset of cognitive and functional decline can be challenging, time-consuming and expensive.
- Episodic memory tasks such as CANTAB Paired Associates Learning (PAL), are associated with AD biomarkers, hippocampal and temporal-frontal network function (de Rover 2011, Nathan 2017) and may therefore be useful for pre-screening.
- Here we assess the utility of PAL to:
 - enrich samples for amyloid positive participants, and
 - identify deficits in other commonly used cognitive tasks such as CERAD learning and delayed recall.

Methods

- Data were evaluated from EDAR, a Europe-wide prospective longitudinal study assessing the potential of biomarkers including cognitive assessments, in the early diagnosis of Alzheimer's disease.
- Participants in the EDAR study comprise a mixed population ranging from normal cognition through MCI to AD and other dementia.
- Estimates of biomarker enrichment were based on the sensitivity and specificity of the CANTAB Paired Associates Learning (PAL) and the age specific probability of Amyloid positivity (Aβ+) in the general population (Jansen, 2015).
- The ability of PAL to detect a deficit in the CERAD was based on logistic regression adjusted by age and gender. A one standard deviation deficit on the CERAD list learning task was calculated using age adjusted standardised norms (Morris, 1989).

Results

- The analysis sample included n=211 participants (Mean age 68.4, SD 9.3) who completed CANTAB and had available biomarker data.
- Participants were divided into Amyloid + (Aβ+) n=169, and Amyloid - (Aβ-) n=42 based on CSF Aβ42 levels (500 pg/ml) assayed using Luminex.

A. Enrichment for Amyloid Positivity

- Figure 1 illustrates the changes in sensitivity and specificity obtained at various thresholds for PAL (dashed lines A-G).
- The figure indicates the proportion correctly classified (true positives (TP), true negatives (TN)) together with misclassifications (false positives (FP) and false negatives (FN)) at each threshold of PAL.

- In this instance, false negatives equate to Aβ+ participants who by not exceeding the number of errors for the PAL cut-off would not be further screened.
- With increasing impairment on PAL, the number of participants required for screening to identify Aβ+ is reduced. The percentage reduction is dependent upon the balance between sensitivity and specificity, together with the pre-test probability. (Figures 2-3)

B. Enrichment for List Learning Tasks

- The correlation between PAL and the CERAD list learning task was -.65 and -.69 for the delayed recall task.
- Performance on PAL was highly specific for identifying a deficit in the CERAD, with an AUC of .89.
- PAL predicted a deficit in CERAD list learning (-1 SD) worse than age adjusted norms with an accuracy of 86% (Figure 4).

Conclusions

- Pre-screening on episodic memory can enrich samples for Aβ+ subjects, and identify those likely to show deficits in other cognitive tasks.
- Scoring algorithms can be tailored according to trial requirements and the protocol specific patient population.
- The CANTAB web-based platform provides a validated and highly sensitive cognitive assessment for pre-screening that can be administered on a home computer and can help reduce the rate of in-clinic screen failures (Cormack, 2016).




Figure 1 - Relationship between AB 42 (y axis) and CANTAB PAL (x axis). A-G relate to thresholds for PAL with increasing impairment from left to right. TN=True Negative, FN=False Negative, FP=False Positive, TP=True Positive.




Illustration (left)
Visual episodic memory assessed using the CANTAB Paired Associates Learning (PAL) task




Figure 2 - Sample Enrichment by PAL. Post-test probability (i.e. enrichment) is calculated from the test sensitivity and specificity together with the pre-test probability (i.e. general population prevalence).




Figure 3 - Illustration of potential reduction in CSF samples required with increasing impairment on PAL. Based on pre-test probability of 29% (Jansen 2015) 345 subjects would be required to obtain 100 positive screens.





Figure 4 - Graph shows the relationship between performance on PAL (y axis) and CERAD list learning (x axis). The vertical line indicates the 1SD cut-off. Coloured dots=predicted group membership logistic regression (red =predicted deficit, blue no deficit) Overall accuracy = 86% Correct, 14% Misclassifications.



Above - Completing CANTAB's Web-based software

References

Cormack F, Taptiklis N, Curtis J et al (2016) Computerised testing: applications in recruitment and monitoring of patients in clinical trials, *J Prev Alzheimer's Dis* 2016;3(4):374-5
 De Rover, M., Pironi, V. a, McCabe, J. et al. (2011) Hippocampal dysfunction in patients with mild cognitive impairment: a functional neuroimaging study of a visuospatial paired associates learning task. *Neuropsychologia*, 49(7), 2060-70.
 Jansen WJ, Ossenkoppele, R, Knol, DL et al. (2015) Prevalence of cerebral amyloid pathology in persons without dementia: a meta-analysis. *JAMA*, doi: 10.1001/jama.2015.4668.
 Morris, JC, Heyman, A., Mohs, R et al (1989). The Consortium to Establish a Registry for Alzheimer's Disease (CERAD). Part 1. Clinical and neuropsychological assessment of Alzheimer's Disease. *Neurology*, 39, 1159-1165
 Nathan P, Lim YJ, Abbott RA et al (2017) Association between CSF biomarkers, hippocampal volume, and cognitive function in patients with amnesic mild cognitive impairment (MCI). *Neurobiology of Aging*, 53

Contact Dr Rosemary Abbott, Email rosemary.abbott@camcog.com; Cambridge Cognition, Tunbridge Court, Tunbridge Lane, Bottisham, Cambridge, CB25 9TU, UK

RAA, FC, KZ, JB, AB are employees of Cambridge Cognition. The authors report no conflict of interest for this work.

EDAR was funded by the European Commission 6th Framework Programme, and the work has been undertaken as part of the European Medical Information Framework (EMIF)



efpia



CAMBRIDGE
COGNITION